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[www.sanowell.com](http://www.sanowell.com)

### **Thank you very much and congratulations**

We are very happy, that you chose a product from hofmann gmbh, and hope, that the sanowell Life Tens helps to relieve your pain.

Please become knowledgeable about your TENS-device before treatment by reading this instruction manual carefully and completely.

The sanowell Life Tens is an analogue TENS device with two controllable output channels for a comfortable therapy in your home. It produces electrical impulses, that may be varied through the control elements in intensity, frequency and width. The entire unit is built very user friendly. Nevertheless, we do recommend prior to therapy to consult your physician.

## 1. General description

The trans-cutic (through the skin) electrical nerve stimulation – short TENS – is a therapy method that has been used extremely successfully for decades to treat intense and chronic pain. With the sanowell Life Tens you have chosen a two-channel-TENS, that will ease or eliminate your pain through electrical impulses, that may be individually set in duration, intensity, frequency, for comfortable treatment at home.

## Indication

Your physician may prescribe a TENS device for you. Please keep strictly to the doctor's orders and have him consult you in correct application of the electrodes. The correct placement is crucial to the success of your TENS therapy.

## 2. What is pain?

Under regular circumstances pain has a warning function within the human body, as it brings our attention to the fact, that something is not okay. Pain is important, as without damage to our body system would not be discovered and therefore not treated. On the other hand pain especially the chronic form puts an immense stress on the body and naturally the human soul. Nobody can and nobody should have to withstand pain for a duration of time.

Pain arises from the brain. Via sensory nerve fibres of our nervous system diversest information reach our brain to be interpreted.

## 3. What is TENS and how does TENS work?

The shortcut TENS (trans-cutic electrical nerve stimulation) depicts the stimulation of nerves by electrical impulse through the skin. In fighting against pain it offers an alternative a-thousand-fold proven in comparison with drugs, that procure a chemical intervention on the body and by that often an additional stress.

The TENS therapy bases on the natural body-own system of information transportation. Positive, as well as negative sensations (as pain) reach the brain by the same nervous system. TENS uses the fact, that comfortable sensations travel on faster nerve strands than negative ones. That

means, when the soft electrical impulses, that TENS-therapy sends out, reach the spinal cord, there is no more room left for the pain information; it can no longer be processed by the brain. Moreover TENS effects that body-own happiness hormones are spilled, the so called endorphins. These two factors together effect that the sensation of pain often is noticeably eased or completely eliminated.

A TENS-therapy may help you to get rid of your pain, so that you can face everyday business more easily. Though do consult your physician and besides pain relief by TENS occupy yourself with exploring the source.

#### **4. What are the fundamentals of treatment?**

##### **4.1. Pulse width**

Pulse width tells how long the duration of one pulse is. Through variation of intensity and pulse width it is possible to stimulate different nerve groups. So it is recommended to use different pulse width settings according to treatment mode and treated area. You may inform yourself about pulse width settings in chapter 5.2.2.

##### **4.2. Frequency**

The frequency of the impulse in Hertz (Hz) tells in what time interval the electrical impulses are sent to the electrodes and thus to the skin. For the choice of frequency the mode of electrode placement is to be considered. In case of stimulation in the pain area or very close to it a higher frequency (80 Hz and up), where the single impulses can not be felt anymore but treatment feels like a continuous tingling, is effective. In case of stimulation of so called trigger points (points in the muscles, that react extremely sensible) or acupuncture points a frequency of 10 Hz or less, where you can feel the single impulse, is recommended.

In combination of both methods higher frequency settings are recommended. Sometimes variations of the recommended settings can be effective. You should consult your physician or try out for the most suitable setting.

You may read about setting the impulse frequency in chapter 5.2.3.

### 4.3. Intensity

Intensity controls the strength of the single impulse. Each patient senses the intensity individually due to tissue resistance, thickness of skin, and so on.

Therefore intensity should be set in a range, where impulses are felt as tingling and tapping, but never as uncomfortable or even painful. Stimulation should not produce muscle contractions. These can be useful when pain is caused by cramped or spastic muscles. Under these circumstances TENS may be used as a common muscle stimulator to release the spasm fast. There it is recommended to set frequency high, pulse width long and intensity to visible contraction (but still tolerable for the patient).

Usually 20 to 30 minutes of such muscle contraction releases the spasm.


Settings for intensity can be found in chapter 5.2.4.

### 4.4. Treatment duration

TENS devices are commonly used for longer periods of time, 20-30 minutes. For some post-surgical treatments the duration may be prolonged to 36 hours. Pain relief may already occur at treatment begin. In some cases it may take up to 30 minutes until pain relieve sets in, especially when point placement of electrodes is applied or at lower frequency settings. Important is to continue pain therapy patiently and consequently.

You may read about settings for treatment duration in chapter 5.2.5.

## 5.0. Draft and arrangement of control elements

### 12.1. Explanation of signs on the CE label

Observe the instruction manual!

Application component of type BF. Application components are not suitable to be used for direct application at the heart.

Symbol for the marking of electric appliances and electronic devices.

## **5.2. Control elements and their function**

### **5.2.1. Setting into operation**

For setting your sanowell Life Tens into operation turn the intensity dial clockwise until you hear a „click“ sound. By turning counter clockwise you can turn off the device.

### **5.2.2. Setting the pulse width**

All pulse width-adjustments are made via the left turning knob (pulse width). The pulse width can be continuously adjusted from 30 to 260  $\mu$ s.

Recommended setting: 200  $\mu$ s

Function of pulse width is explained in chapter 4.1.

### **5.2.3. Setting impulse frequency**

All pulse frequency-adjustments are made via the right turning knob (pulse rate). The pulse frequency can be continuously adjusted from 2 to 150 Hz.

Recommended setting: 100 Hz

Function of impulse frequency can be read in chapter 4.2.

### **5.2.4. Setting intensity**

Intensity is set by turning the intensity dials. Turning clockwise increases, turning counter clockwise decreases intensity steadily. Always start with the lowest intensity setting and work your way up slowly until you feel a tingling sensation. Treatment never ought to be uncomfortable or even painful.

The function of intensity is described in chapter 4.3.

### **5.2.5. Setting treatment duration**

With the help of the lower slide control (timer min) you can choose between three treatment durations: C for the continuous mode, 30 or 60 minutes.

Recommended treatment duration: 30 minutes

The function of treatment duration can be read in chapter 4.4.

## 6. Preset programs

With the help of the upper slide control (mode) you can choose between three preset programs: B (=burst), C (=continuous) or M (=modulation).

Setting suggestion: „M“

### 6.1. B, or burst treatment method

In the first treatment method, burst (B), single bursts of 7-10 individual impulses in a fixed pattern are sent. This is a combination of common TENS and low frequency TENS. In the burst mode the frequency of stimulation is set by the device and can be controlled separately through the frequency control.

### 6.2. C or continuous treatment method

The medium adjustment for a continuous treatment method (C), allows the user to diversify the intensity, pulse width and pulse rate manually.

How to conduct setting can be found in chapter 5.

### 6.3. M, or modulated treatment method

The third setting is the modulated treatment mode (M). through continuous change in intensity this method tries to prevent immunisation of the nerve. If intensity is increased in a period of low intensity within the cycle, intensity can be increased slowly by the control element, so that you will feel this as an increase.

## 7. What to observe about the electrodes

Your physician ought to decide, which electrodes to best use for your kind of pain. Do observe the instructions on electrode placement in the instruction manual in order to keep stimulation effective during treatment time and to prevent skin irritations. You may only use electrodes with the CE-certification sign.

The wires delivered with the TENS device are plugged into the sockets at the device. To do so place the device in your palm and push the plugs at the end of the wire into one of the sockets (see drawing); you may use one or two connection wires. After connecting the wires please connect the electrodes to the wire ends (positive and negative terminal). Please be careful when attaching or disconnecting the wires. Do not pull at the wires, as this may lead to breakage.

## **7.1. Electrode placement**

Best possible electrode placement is the key to successful treatment. Besides the most usual placement strategies described as follows, it may be effective to use a deflective electrode application. Thus it is important that you take your time in experimenting for the best electrode placement for you, as the positioning is individual. Please have a TENS firm physician consult you.

### **7.1.1. The bordering placement**

This is the most common placement strategy. Here the painful area is bordered by four electrodes. You can decide whether the current is to flow directly through the pain area or surrounds it. In order to obtain a direct flow through the spot of pain you have to place the positive terminals of the electrode pairs diagonally to the negative ones; if trying to achieve a surrounding flow, they have to be placed parallel. This is very effective if the pain is in an extremity or deep within the tissue.

### **7.1.2. Placement in the area of a spinal nerve**

The spinal nerves exit the spinal cord in pairs between the vertebrae. In this technique of electrode placement stimulation happens above the innervated area by placing one electrode exactly on the spot of pain and the other at the nerve's root at the spinal cord.

### **7.1.3. Placement at motive, trigger, and acupuncture points**

Motive, trigger, and acupuncture points are places of high tissue conductivity. In order to destine the exact position of such a point please consult your physician.

### **7.1.4. Multiple placement strategies**

As your TENS device possesses two independent channels, you can apply different strategies in electrode placement at the same time. One channel can be used for direct stimulation of the pain area, whereas the second channel is applied for point therapy by bordering placement.

## **7.2. Use of self-adhesive electrodes**

### **7.2.1. Application**

- Clean the skin thoroughly with soap and water before application of the electrodes

- Excessive hair should be cut, but not shaved
- Connect the wires to the electrodes
- Pull the safety cover off the electrodes and apply firmly to the skin area by pressing.

### **7.2.2. Removal**

- Lift the electrodes at the corner and gently pull off. Do not pull at the electrode wire, as this might damage the electrodes
- Re-glue the electrodes to the safety cover and then remove the wires by carefully turning and pulling.

### **7.2.3. Care and storage**

Your sanowell Life Tens will only achieve best results when electrodes and wires are cared for accordingly.

- Clean the wires with a damp cloth; slight powdering with talcum (f. ex. baby powder) will prevent tangling of the wires and improves their durability
- After each application please store the electrodes in their re-closeable bags in a cool dry place (f. ex. the refrigerator)
- Life time of the electrodes is improved by applying a drop of cool tap water (please do not use distilled water) to the adhesive area and air drying the surface turned up. Moistening the adhesive area too much may spoil the adhesive capability.

## **8. General safety and precautionary measures**

Do read the instruction manual, before using your sanowell Life Tens for the first time.

- Only use your TENS device according to its purpose described in the instruction manual.
- Only employ accessories recommended by hofmann gmbh.
- Persons with electronic or metallic implants (f. ex. pace maker) may not be treated by TENS without prior permission of their physician.
- In case TENS treatment is ineffective or uncomfortable stimulation has to be ceased until examination by a physician has been conducted.
- Do not apply the device when operating machinery or driving a vehicle.
- Please always turn off your device before application and removal of the electrodes.

- TENS devices do not have a AP/APG-safety. They may not be operated in the environment of explosive or combustible material.
- Persons afflicted with heart disease may only apply the sanowell Life Tens under extreme caution. Please consult your physician concerning side effects before application of the device.
- TENS may not be applied during pregnancy.
- The purposeful use of TENS is nerve stimulation. If the device is applied to unsuitable body areas danger to the user may arise. Never place electrodes in the laryngeal or pharyngeal area. Muscle spasms in this area can afflict asphyxiation.
- In some rare occasions skin irritations may be caused by the electrodes. In such incidence please cease stimulation and remove the electrodes. Discontinue treatment until the cause for the skin irritation has been diagnosed.
- Do not apply electrodes to the eyes, in the mouth, or in other body holes. The electrodes are to be used only on healthy skin. Never apply electrodes to irritated or injured skin.
- Keep away from children.
- Never apply electrodes in a manner that will trans-duct currency through the brain.
- The device may not be used to treat pain of uncertain cause.
- Attention: Do not plug the electrode wires into a alternating current socket!

### **Additionally in Germany it is applicable as follows:**

- The device is manufactured according to the medical products law (MPG).
- The device is subject to the medical products operator's regulation (MPBetreibV).
- Operator and user are obligated to strictly observe the operator's regulation in all aspects.
- For all other states of the European Community national laws and regulations for operation of medical products are valid.

### **9. Helpful hints for skin care**

In order to prevent skin irritations we suggest as follows:

- Wash the areas where the electrodes are to be placed prior to application and after removal with a mild soap and water. Rinse the soap thoroughly and dry the skin.

- Excessive hair growth may be cut away with the scissors; do not shave the area that is to be stimulated.
- Many skin problems arise due to the electrodes being pulled too tight over the skin at application. To prevent this apply the electrodes from their centre to the rim and do not pull them over the skin.
- To prevent irritation by stretching tape the excess wires in a circle to the skin so it will not pull on the electrodes.
- When removing the electrodes always do so in the direction of hair growth.
- It may be helpful to apply skin lotion on the areas that electrodes were attached to between applications.
- Never adhere electrodes to skin already irritated or injured.

### **Attention**

- Do not apply to injured skin.
- Dispose of electrodes that do not adhere well.
- Self-adhesive electrodes are to be used by the same person only.
- In case of skin irritation cease usage and consult your physician.
- Prior to application read the instruction manual for self-adhesive electrodes.

If stimulation is uncomfortable either reduce intensity to a tolerable amount or discontinue stimulation and ask your physician.

## **10. Handling the Life Tens**

Your sanowell Life Tens is already equipped with a battery, but may also be operated by storage batteries.

### **10.1. Usage of re-chargeable batteries**

Prior to usage of a new device the storage battery has to be charged according to the manufacturer's prescription. Please read all instructions and safety precautions printed on the accumulator and in this leaflet before usage. Stored for more than 60 days the re-chargeable battery may discharge. Therefore prior to usage they need to be charged again, if they were stored for a duration of time.

### 10.1.1. Charging the accumulator

- Plug the accumulator into a 220/230V socket. Using an accumulator not delivered as accessory may cause fire, electric shock or body injury.
- Regarding charging time observe the instructions of the manufacturer.
- After elapse of charging time unplug the accumulator and remove the charged battery.
- Re-chargeable batteries should always be stored loaded.

To guarantee best battery performance the following rules should be observed:

- Overloading of up to 24 hours does not damage the re-chargeable battery. Though frequent overloading may reduce their life time.
- Do not short plug the accumulator terminals. This causes the accumulator to run hot and may cause permanent damage. Do not keep batteries in clothing or purses as terminals may contact coins, keys, or other metal objects.

### Safety precautions:

- Do not charge any other batteries with the accumulator than re-chargeable batteries. Other tape batteries may leak or explode.
- Do not burn re-chargeable batteries, as they may explode!

### 10.2. Testing and exchanging batteries

A 9V battery is already included in delivery. To ensure the device's operational reliability the battery needs to be changed from time to time.

When replacing the battery proceed as follows:

- Pull off the battery compartment's lid downwards.
- Remove the consumed battery from the battery compartment.
- Insert a new 9V battery according to the drawing in the bottom of the battery compartment. The destined polarity has to match.
- Push the lid back upwards.

### **10.3. Care, transportation, and storage**

- The device may be cleaned with alcohol. Attention: When handling combustible fluids do not smoke or operate with open flame (e.g. candles, etc.)
- Stains and massive soiling may be removed with a mild detergent.
- Do not immerse the device into liquids or subject it to large amounts of water.
- For safe transportation always insert the device into the carrying case.
- If the device is not used for a prolonged period of time, remove the batteries from the battery compartment (acid could leak and damage the device). Store the device and all accessories cool and dry in the carrying case.

### **10.4. Safety-technical control measures, STC**

According to the medical products operator's regulation (MPBetreibV), § 6 (safety-technical control measures), the user is obliged to have regular safety-technical controls conducted. hofmann gmbh prescribes these controls in accordance to MPBetreibV § 6 in a 24-months turn.

The safety-technical control may only be delegated to persons, who are able to conduct controls in an orderly fashion based on their education, knowledge, and experience achieved by practical work, and who also have suitable measuring and proofing tools available.

Standards for safety-technical control are to be requested from the quality management agent of hofmann gmbh.

Prior to each application of the sanowell Life Tens device the following controls are to be conducted:

1. Inspect device for surface damage:
  - Deformation of casing?
  - Damaged or defect electrode sockets?
2. Inspect device for damaged control elements:
  - Imprints still legible?

3. Inspect display elements:
  - When the corresponding channel is working the related LED needs to be illuminated.
4. Inspect accessories for usability:
  - Electrode wires non-damaged?
  - Electrodes non-damaged (no fissures, no worn material)?

### 10.5. Malfunctions

Should malfunctions occur, when using the TENS, please check as follows:

- Is the treatment mode set correctly? If not, please correct setting.
- Is the wire connected correctly to the device? Wires are to be completely inserted into the sockets.
- Do both LED bulbs work? Change battery, if necessary.
- Is the wire damaged? In any case of damage, exchange the wire.
- Are 2 Electrodes connected to the paired wire?

In case of any other problems, return the device to the dealer. Do not try to repair a defect device yourself! Manipulation voids the warranty!

For problems with the device or accessories please contact your dealer.

### 10.6. Disposal

Device packaging should be disposed of in the recycling bin, metal parts of the device in the metal recycling disposal. Plastic parts, electrical parts, and electronic chips are to be disposed of as electronic garbage. Batteries are problem garbage and are to be disposed of accordingly. Return consumed batteries to your local recycling plant. Do not dispose of in regular garbage. You may obtain further information on disposal of problem material at your local waste management centre.

## **11. Contents of delivery and accessories**

### **11.1. Contents of delivery**

- 1 x sanowell Life Tens unit
- 2 x lead wires
- 4 x self-adhesive electrodes, size 40 x 40 mm
- 1 x battery, 9 V
- 1 x instruction manual
- 1 x carrying case

#### **Caution:**

The above mentioned contents are standard equipment, that can vary within distributors.

### **11.2. Accessories**

- self-adhesive electrodes, 2mm plug, 40 x 40 mm, article No. 50410
- self-adhesive electrodes, 2 mm plug, 40 x 80 mm, article No. 50100
- self-adhesive electrodes, 2 mm plug, size 32 mm round, article No. 50050
- permanent electrodes for use with electrode adhesive gel, article No. 50430
- electrode gel for TENS and EMS devices, article No. 14221
- cables for permanent electrodes, article No. 50420
- cables for self-adhesive electrodes, article No. 50400

**If you would like to make your device compatible with the permanent electrodes, you need to supply yourself with the fitting cables and the electrode gel.**

## 12. Technical specifications

Channels:	dual, isolated between channels
Pulse Amplitude:	adjustable to maximal 80 mA for positive impulse or to about 10 mA peak for negative impulse at 500 Ohm load per channel.
Pulse Rate:	adjustable from 2 to 150 Hz
Pulse Width:	adjustable from 30 to 260 $\mu$ s
Timer:	30 and 60 minutes or continuous duration
Wave Form:	asymmetrical bi-phasic square pulse
Modulation mode:	Within 6,5 seconds the pulse duration varies cyclic. The duration diminishes linearly within 0,5 seconds to a value that is 60% lower than the initial value. This duration is continuous for 2 seconds afterwards it is increased again within 0,5 seconds up to the initial value which is held for 3,5 seconds before the cycle starts again.
Burst mode:	Bursts consist of pulse width (adjustable) and a frequency of 100 Hz. Bursts are effected twice a second.
Power Source:	9V battery
Battery Life:	approximately 60 hours at regular settings
Operating Conditions:	0 °C to +50 °C at 20 to 75% relative humidity, however with no condensation, pressure 700 hPa to 1060 hPa.
Storage Conditions:	0 °C to +50 °C at 20 to 75% relative humidity, however with no condensation, pressure 700 hPa to 1060 hPa.
Dimensions:	95 mm (height) x 65 mm (width) x 23,5 mm (length)
Weight:	115 g (inclusive of battery)
Safety for electrical shock, type BF.	
All electrical specifications valid at $\pm$ 20 % 500 Ohm load.	

### **13. Warranty**

There is a 2 year guarantee for our products, applicable from date of purchase.

Within the prescribed warranty time we will remove all defects of the device, that are based on faulty material or manufacture. Removal of defect occurs at our choice by either repair or replacement of parts or the entire device.

Exempt from warranty are: Defects caused by inappropriate use (f. ex. connection to inappropriate electrical current or voltage, connection to inappropriate sources of electricity, breakage by fall and such). Normal wear and tear that only minimally affects the operational utility of the unit is not subject to the warranty.

Warranty subsides if unauthorized intervention or repair is undertaken, or if non-original parts are used.

In case of a defect that is covered by the warranty, please send the device and the original receipt to one of our service partners or directly to hofmann gmbh.

No other claims are valid under the provisions of this warranty, unless our liability is based on legal statute. Claims against the dealer based on the provision of the purchase agreement do not affect this warranty.

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Hofmann GmbH