# visomat®

handy

Die  
mobile  
Messung

Gebrauchsanweisung  
Instructions for use  
Mode d'emploi  
Istruzioni per l'uso

DE-2 - DE-29  
EN-30 - EN-55  
FR-56 - FR-83  
IT-84 - IT-111

**UEBE**  
Germany est.1890

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## Table of contents

|          |  |    |
|----------|--|----|
| <b>A</b> | <b>Intended use</b>                                | 32 |
| <b>B</b> | <b>Safety instructions</b>                         |    |
|          | 1. Important patient information                   | 33 |
|          | 2. Important technical information                 | 34 |
| <b>C</b> | <b>Operating the unit</b>                          |    |
|          | 1. Unit description                                | 37 |
|          | 2. Display   | 38 |
|          | 3. Important instructions for use                  | 38 |
|          | 4. Initial operation of the device                 | 39 |
|          | 5. Inserting/replacing batteries, battery icon     | 40 |
|          | 6. Activating and setting the date/time            | 40 |
|          | 7. Attaching the cuff                              | 41 |
|          | 8. Body posture when taking measurements           | 42 |
|          | 9. Measuring blood pressure                        | 43 |
|          | 10. Classification of the measured values (WHO)    | 43 |
|          | 11. Pulse pressure                                 | 44 |
|          | 12. Irregular pulse waves and cardiac arrhythmia   | 45 |
|          | 13. Using the memory                               | 45 |
| <b>D</b> | <b>What you should know about blood pressure</b>   |    |
|          | 1. Systolic and diastolic blood pressure values    | 47 |
|          | 2. Reasons for measuring different values          | 47 |
|          | 3. Why you should measure blood pressure regularly | 47 |
| <b>E</b> | <b>Technical information</b>                       |    |
|          | 1. Error messages                                  | 49 |
|          | 2. Customer service                                | 50 |
|          | 3. Technical data                                  | 51 |
|          | 4. Original spare parts and accessories            | 52 |
|          | 5. Metrological inspection                         | 52 |
|          | 6. Explanation of symbols                          | 53 |
|          | 7. Disposal  | 53 |

## Table of contents

|          |                             |    |
|----------|-----------------------------|----|
| <b>F</b> | <b>Maintaining the unit</b> | 54 |
| <b>G</b> | <b>Warranty</b>             | 55 |

Thank you for choosing the visomat® handy wrist blood pressure monitor (also referred to in the following as the unit).

These instructions for use are intended to help the user operate the device safely and efficiently. The unit must be used in accordance with the procedures contained in these instructions for use and must not be used for any other purposes.

It is important that you read all the instructions carefully before using the device. Please pay particular attention to the chapter "Important usage information" on page 38.

**A****Intended use**

This device is intended for persons aged 12 or over for non-invasive self-measurement of the systolic and diastolic blood pressure, determining the pulse rate and calculating the pulse pressure.

This unit uses the oscillometric method for measuring blood pressure and pulse rate.

**Safety instructions****B****1. Important patient information**

- The device may only be used at the wrist and not on any other extremities.
- Blood pressure measurements on children require specialist knowledge! Please consult your doctor if you wish to measure the blood pressure of a child. Do not under any circumstances use the unit on a baby/infant.
- During pregnancy and in cases of pre-eclampsia, the device may only be used after consulting your doctor.
- Do not under any circumstances place the cuff over a critical area, e.g. a wound, aneurysm, etc. or on an arm with an arteriovenous shunt. Risk of injury! Any supply via an intravascular access point (infusion) or other medical monitoring devices could possibly be interrupted.
- Do not use the unit without first consulting your doctor if you are undergoing dialysis treatment or taking anticoagulants, platelet aggregation inhibitors or steroids. Internal bleeding may be caused under these conditions.
- The display of the pulse frequency is not suitable for checking the frequency of cardiac pacemakers. Cardiac pacemakers and blood pressure monitors do not influence each other regarding their mode of operation.
- The device is not suitable for use in close proximity to electro-surgical instruments or magnetic resonance imaging equipment.
- Do not use the device in explosive environments such as in the vicinity of flammable anaesthetics or in an oxygen chamber.

- Please wait a few minutes between successive measurements, otherwise the blood flow in the arm is interrupted for too long and injuries could occur.
- The device contains small parts that could be swallowed by children. Do not leave the device unattended with children.
- Store the device out of the reach of children.
- Do not take a blood pressure measurement on the same side as a mastectomy with simultaneous removal of the axilla lymph nodes.
- Discuss the measured values with your doctor. Do not assess the measurement results yourself. **Never change the dose of medicines prescribed by your doctor.**
- Before carrying out your self-measurement, please pay attention to the chapter "Important usage information" on page 38.

## 2. Important technical information

- The unit contains sensitive parts and must be protected from extreme temperature fluctuations, humidity, shocks, dust and direct sunlight.
- Do not knock or drop the device.
- The device is not waterproof. Never immerse it in water or any other liquids. If liquid enters the device, this may result in serious damage and may cause malfunctions.
- Do not apply pressure to the display. Do not place the device with the display facing downwards.
- The unit must only be operated with the cuff and accessories provided (see "Original spare parts and accessories" on page

- 52). Using different cuffs or accessories can lead to incorrect measurement results. The warranty will be invalidated if the unit is damaged by third-party accessories!
- A consistently good power supply to the device is necessary for fault-free blood pressure measurements.
  - Only use long-life alkaline batteries (LR03).
  - Always replace all batteries simultaneously during battery replacement. Do not mix new and old batteries or batteries of different types.
  - You will require 2 x 1.5 V batteries. Rechargeable batteries only have a voltage of 1.2 V and are therefore not suitable.
- Never simultaneously touch the device's live electrical parts and the patient.
- If the device is not to be operated for an extended period of time, please remove the batteries. In principle, all batteries can leak.
- Do not remove batteries from the device when the device is switched on.
- Using this device in the vicinity of mobile phones, microwaves or other devices with strong electromagnetic fields can cause malfunctions and inaccurate measurements. During use, maintain a distance of at least 30 cm from such devices.
- Never open or modify the device or the cuff (except for replacing the batteries). If the unit has been opened, it must be subjected to a metrological inspection by an authorised institution.
- To avoid inaccurate readings, please adhere to the intended operating and storage conditions. See "Technical data" on page 51.


- The inflating and measuring procedure can be interrupted by pressing the Start/Stop button or by removing the cuff. The device then stops the inflation procedure and deflates the cuff.
- If the device malfunctions or is faulty, please refer to the troubleshooting guide starting on page 49, or contact customer service (see page 50).

## 1. Unit description

- |                     |   |
|---------------------|---|
| 1 Display           | 4 Memory button   |
| 2 Cuff              | 5 Battery compartment (Removable battery compartment cover) |
| 3 Start/Stop button |   |

2. Display

- 1 SYS = Systole
- 2 DIA = Diastole
- 3 PUL 1/min = Pulse Calculated pulse frequency per minute
- 4 Pulse signal display or irregular pulse waves
- 5 Battery check display
- 6 Inflating / Deflating
- 7 WHO classification (Page 43)
- 8 Memory identifier
- 9 Storage location or date/time

3. Important instructions for use

The measurement results of automated blood pressure monitors can be influenced by the measuring location, posture, previous exertion and general physical condition. Please pay attention to the usage information in order to obtain accurate readings.

- Do not consume alcohol, caffeine or smoke for at least one hour before measuring.

- Rest for at least 5 minutes before taking the measurement. Depending on the previous degree of exertion, you may even need to rest for up to an hour.
- Expose the wrist. Under no circumstances may clothes hinder the blood flow to or from the arm, as this affects the blood pressure at the measuring point and can lead to incorrect readings.
- Do not move or talk during measuring.
- Breath slowly and deeply. Do not hold your breath.
- Check the display for irregular pulse signals. If required, repeat the measurement when conditions have improved.
- In order to exclude side differences and to obtain comparable measurement results, it is important to always take the measurement on the same arm. Ask your doctor which side it is advisable for you to measure.
- Always take your daily blood pressure measurements at around the same time.
- Blood pressure is not a fixed value. It can go up or down by more than 20 mmHg in patients within a few minutes.

4. Initial operation of the device

Insert the batteries provided into the device.

The device has a date/time function that can be activated as required (see page 40). This function is deactivated on delivery.

When the date/time function is activated, the measurement results are saved along with their date and time. This can make it easier to view the readings at a later point.

## 5. Inserting/replacing batteries, battery icon

- Remove the cover of the battery compartment on the side of the unit.
- Remove the old batteries from the device and insert the new ones. Please ensure correct polarity (marking in battery compartment).
- Close the battery compartment by clipping the battery cover back into the device.

When the “empty battery ” symbol first appears in the display, you are still able to take around 30 measurements. Please replace the batteries within this time frame.

## 6. Activating and setting the date/time

### Activating and deactivating the date/time function

To activate the date/time function, press the Start/Stop button and then immediately, while the full display is visible, the Memory button. This function can be deactivated at any time in the same way.

When the date/time function is activated, the device continues to display the time even when it is switched off. The power consumption for this is extremely low.

### Setting the date/time

The four-digit number representing the year appears in the display first (Figure 1). Press the Memory button to set the year. Press the Start/Stop button to confirm the selected year and proceed to setting the date (Figure 2).

Figure 1

Set the month with the Memory button and confirm this with the Start/Stop button. Set the day, hours and minutes in the same way (Figure 2, Figure 3).


Figure 2


Figure 3

### Changing the date/time


To change the date and time, take the batteries out and wait until the display goes off. Then re-insert the batteries and reset the date and time.

After replacing the batteries, you will have to reset the date and time.

## 7. Attaching the cuff

The special double-preformed comfort cuff is easy to attach and ensures a perfect fit during measurement. The materials used are tested for skin compatibility according to ISO 10993.

- Clear the wrist from all clothing or watches/jewellery since this may disturb the blood flow in the arm as it may impair the blood pressure at the place of measurement.
- Slide the device onto the centre of the inside of the forearm, directly on the wrist. The palm of your hand should be facing upwards. The distance between the wrist and the upper edge of the cuff should be 5-10 mm.
- Fasten the cuff using the hook and loop fastening.


The shape of the cuff is optimized to sit very tightly in order to avoid possible measurement errors.

## 8. Body posture when taking measurements

The measurements should be performed while seated. In special cases, measurements in a lying position are also possible.


### Measurement while sitting

Sit at a dining table, lean your back against the backrest of the chair and bend your arm so that your arm is exactly at heart level. Place your feet on the ground and do not cross your legs.

### Measurement while lying down

Lay on your back, place the device on your wrist and at the height of your heart, e.g. by supporting it with a pillow. To relax before the measurement, take 2-3 deep breaths.

### Please pay attention to the following:

- Do not bend your wrist.
- Keep your hand relaxed, not tightened, do not make a fist.
- It is important to keep your arm calm and not to speak during the measurement.

If the cuff is in a higher or lower position than the heart, the measured value might be higher or lower, respectively.

The results for measurements taken in a sitting and lying position may differ slightly.

## 9. Measuring blood pressure

- Switch on the unit by pressing the Start/Stop button. The full display appears (Figure 1).
- The cuff is inflated to the required pressure. The ▲ symbol flashes in the display (Figure 2).
- The actual measurement process starts with three beeps as the cuff is deflated. The ♥ symbol flashes in the display (Figure 3).
- A long beep denotes the end of measurement. The cuff is automatically deflated.
- The determined values of systole, diastole and pulse and the pulse pressure (PP) are shown alternately in the display (Figure 4, Figure 5).
- The device switches itself off after approximately 3 minutes. You can also switch the device off before this by pressing the Start/Stop button.


Figure 1

Figure 2

Figure 3

Figure 4

Figure 5

## 10. Classification of the measured values (WHO)

- Upper dot (Red) signifies high pressure\*
- Middle dot (Yellow) signifies borderline value\*
- Lower dot (Green) signifies normal value

| WHO 2003 | Systolic pressure<br>= Upper value<br>mmHg** | Diastolic pressure<br>= Lower value<br>mmHg** |
|----------|--|---|
| Red      | from 140*                                    | from 90*                                      |
| Yellow   | 120 to 139                                   | 80 to 89                                      |
| Green    | below 120                                    | below 80                                      |

\* It is sufficient for one of the values to be raised.

\*\* millimetres of mercury

The categorisation according to WHO is stored with the readings and can be called up again together with the readings from the memory.

Depending on age, weight and general health, the blood pressure values can be different. Only a doctor can determine the right blood pressure range for you and assess whether your blood pressure has reached a dangerous level for you. Discuss your blood pressure values with your doctor. **Never change the dose of medicines prescribed by your doctor!**

## 11. Pulse pressure

The pulse pressure – not to be confused with the pulse beat – provides an indication of the elasticity of the blood vessels. A stiff system of vessels can have a negative effect on the cardiovascular system. Studies show that the cardiovascular risk increases if the pulse pressure value is permanently higher than 65 mmHg.

The heart works in two phases, the contraction phase (systole) and the relaxation phase (diastole). The pressure difference between systole and diastole is called the pulse pressure or pulse amplitude. The higher the pulse pressure, the harder the vessels.

|                          |               |
|--------------------------|---------------|
| High pulse pressure      | above 65 mmHg |
| Increased pulse pressure | 55 to 65 mmHg |
| Normal pulse pressure    | below 55 mmHg |

If your pulse pressure is permanently above 65 mmHg, consult your doctor.

## 12. Irregular pulse waves and cardiac arrhythmia

If the device detects irregular pulses during the measurement, the  symbol is displayed once the measurement is complete. This can be set off by an irregular heartbeat (arrhythmia), disturbances caused by movement, talking or even breathing in deeply. The symbol is stored with the respective measurement.

If the symbol appears more frequently, this may be caused by an irregular heartbeat and must be discussed with your doctor!

Measurement results accompanied by irregular pulses should therefore be viewed as critical and repeated under more favourable conditions.

**Caution:** Do not confuse the irregular pulse display  with the flashing pulse display  that indicates the heart rate during every measurement (see page 43).

## 13. Using the memory

Measured results are automatically stored in the memory. The memory can store up to 60 results and the average value.

When more than 60 measured values have been stored, the oldest value (No. 60) is deleted to allow the latest value (No. 1) to be recorded.

### Data retrieval

Press the Memory button to call up data. The average value of the stored results is displayed with "A", the values for systole, diastole and pulse and the pulse pressure (PP) are shown alternately in the display (Figure 1, Figure 2).


Figure 1 Figure 2

When the Memory button is pressed again, the value measured last (No. 1) appears. When the date/time function is activated, systole, diastole and pulse are displayed alternately with memory number (Figure 3) and date (Figure 4), followed by the pulse pressure (PP) and the time (Figure 5).


Figure 3 Figure 4 Figure 5

Press the Memory button repeatedly to call up further measured values.

The stored data remain visible in the display for around 30 seconds. The unit then switches itself off.

### Deleting data

To delete individual measured values, press the Memory button repeatedly to call up the desired measured value.

Then press and hold the Memory button again for 8–10 seconds until the measured value flashes several times and then goes out.

To delete the entire memory, retrieve the average value then press and hold the Memory button again for 8–10 seconds until the measured value flashes several times and then goes out.

If no value is stored, there is no display.

## 1. Systolic and diastolic blood pressure values

The cardiovascular system has the important function of supplying all organs and tissues in the body with sufficient amounts of blood and of transporting metabolites. For this, the heart contracts and expands at a regular rate of about 60 to 80 times per minute. The pressure of the flowing blood on the artery walls caused by the heart contracting is termed systolic. The pressure in the ensuing relaxation phase, when the heart refills with blood, is termed diastolic. During daily measurement you determine both values.

## 2. Reasons for measuring different values

Our blood pressure responds to internal and external influences like a sensitive measuring instrument. Affected by mental, physical and environmental influences, it varies all the time and never remains constant.

Reasons for fluctuating blood pressure values can be: moving, speaking, eating, consuming alcohol or nicotine, nervousness, inner tension, emotions, room temperature, recent urination or bowel movements, environmental influences such as movements and sounds, etc. Even changes in the weather and climate can have an effect on your blood pressure.

This explains why values measured at the doctor are often higher than those you obtain at home in your usual environment.

## 3. Why you should measure blood pressure regularly

Even the time of day has an effect on your blood pressure. The values are at their highest during the day. In the course of the afternoon and in the evening, they drop slightly. They are low while you are sleeping, but rise again relatively quickly once you get up.

## D What you should know about blood pressure

One-off and irregular measurements therefore say little about your actual blood pressure. A reliable assessment is possible only when measurements are taken regularly. Discuss the measurement values with your doctor.

## Technical information

### 1. Error messages

| Failure encountered   | Possible cause  | Corrective action   |
|---|---|---|
| Display<br>Err - 300  | Excessive cuff pressure. Cuff inflated to maximum due to movement of the arm or body.             | - Repeat the measurement<br>- Do not move your arm<br>- Do not speak  |
| Display<br>Err - 1  | Measurement error; measurement could not be carried out due to movement or speaking.              | - Repeat the measurement<br>- Do not move your arm<br>- Do not speak  |
| Display<br>Err - 2  | Inflation error, pressure built up too quickly or too slowly.                                     | - Repeat the measurement<br>- Do not move your arm<br>- Do not speak  |
|   | Cuff is not correctly fitted.   | - Check wrist circumference. Call the Service hotline if this occurs repeatedly.  |
| Display<br>Err - 3  | Deflation rate too high or too low.   | - Repeat the measurement<br>- Do not move your arm<br>- Do not speak  |
|   | Movement or talking during measurement.   | Call the Service hotline if this occurs repeatedly.   |
| Display 0<br>0  | Start/Stop button accidentally pressed while batteries were being changed.                        | Switch unit off and then on again with Start/Stop button.   |
|  | Irregular heartbeat, disruptive movements, shaking, wobbling (objects), breathing in deeply, etc. | Repeat measurement after 3-5 minutes rest. Possible influences on the measurement results (depending on the severity of the arrhythmia) must be discussed with your doctor! |
| Measured values/<br>pulse values are<br>extremely high<br>or low.                 | Too little rest before the measurement. Movement or talking during measurement.                   | Repeat the measurement after a break of 3-5 minutes. Pay attention to the usage information on page 38.   |
|   | Incorrect measurement posture   | Repeat the measurement and maintain the correct posture. (Page 42)  |

| Failure encountered  | Possible cause   | Corrective action   |
|--|--|---|
| Consecutive measurements produce different measured values (readings).   | Blood pressure is not a fixed quantity. It may vary in the upward or downward direction within a few minutes by more than 20 mmHg. | In order to obtain comparable results, always carry out the measurements under the same conditions. |
| Measurement was interrupted, the cuff deflated to 0 and the measurement restarted.   | Device detects a fault or that the diastolic value is too low.   | Device shows no fault. Repeat measurement   |
|  | Movement during measurement  | If necessary, interrupt measurement, measure again after a break of 5 minutes.                      |
| No display after unit is switched on   | Batteries inserted incorrectly.  | Check position of batteries.  |
|  | Batteries flat.  | Change batteries.   |
|  | Battery contacts dirty.  | Clean battery contacts with a dry cloth.  |
| Nothing appears in the display after pressing the Memory button  | No measured values have been stored.   | Repeat measurement  |
| Measurement interrupted  | Batteries flat.  | Change batteries.   |
| The  symbol appears during measurement but the measurement continues. | Low batteries  | You can still perform approx. 30 measurements, but should keep spare batteries available.           |
| Date/time not shown in display   | Function accidentally deactivated  | Activate date and time as described on page 40.   |

## 2. Customer service

Device repairs may only be performed by the manufacturer or an expressly authorised body. Please contact:

UEBE Medical GmbH  
 Zum Ottersberg 9, 97877 Wertheim, Germany  
 Tel.: +49 (0) 9342/924040, Fax.: +49 (0) 9342/924080  
 Email: info@uebe.com, Website: www.uebe.com

## 3. Technical data

|  |   |
|--|---|
| Model:                                   | visomat® handy  |
| Size:                                    | approx. W = 70 mm x H = 70 mm x D = 80 mm with cuff   |
| Weight:                                  | 121 g with cuff, not including batteries  |
| Display:                                 | LCD display (liquid crystal display) for measured values and check displays   |
| Memory:                                  | 60 measured values (stored automatically) and average value (A)   |
| Measurement procedure:                   | Oscillometric determination of systole, diastole and pulse  |
| Reference procedure of clinical testing: | Auscultatory measurement  |
| Inflating pressure:                      | Automatically controlled  |
| Pressure display range:                  | 0-300 mmHg  |
| Measurement range:                       | Systolic: 50-250 mmHg<br>Diastolic: 40-180 mmHg<br>Pulse measurement: 40-160 pulses/min                             |
| Measurement precision:                   | Blood pressure measurement (cuff pressure):<br>± 3 mmHg<br>Pulse rate: ± 5%   |
| Power supply:                            | Battery type: 2 x 1.5 V, AAA, LR 03, durability: over 700 measurements in 2 years                                   |
| Cuff:                                    | Cuff type visomat H5, wrist circumference 12.5 cm - 21.5 cm   |
| Operating conditions:                    | Ambient temperature 10°C to 40°C, relative humidity 15% to 85% non-condensing, air pressure 700 to 1060 hPa         |
| Storage and transportation conditions:   | Ambient temperature -5 °C to 50 °C<br>Relative humidity up to 85%, non-condensing                                   |
| Inflation and deflation:                 | Regulated automatically   |
| Protection class IP:                     | IP 20: Protected against solid foreign particles with a diameter of more than 12.5 mm, no protection against water. |
| Protection against electric shock:       | Internal power supply, applied part type BF (cuff)  |
| Expected service life:                   | 5 years   |

|                 |                                 |
|-----------------|---------------------------------|
| Operating mode: | Continuous operation            |
| Classification: | Internally powered ME Equipment |

#### 4. Original spare parts and accessories

The following original spare parts and accessories are available from specialist dealers:

- Cuff type H5 (for wrist circumference 12.5-21.5 cm)  
Part no. 2105001, PZN-06414493

Subject to technical modifications.

#### 5. Metrological inspection

Generally, it is recommended that a metrological inspection is performed every two years. However, professional users in Germany are obligated to do so in accordance with the "Medical Products Operator Ordinance" (Medizinprodukte-Betreiberverordnung). This can be performed either by UEBE Medical GmbH, an authority responsible for metrology or an authorised maintenance service. Please refer to your national regulations.

Upon request, responsible authorities or authorised maintenance services receive a "Test instruction for metrological inspection" from the manufacturer.

Please only submit the device for metrological inspection together with the cuff and the instructions for use.

Important: No modifications, e.g. opening the device (except to replace the batteries), may be made to this device without the manufacturer's permission.

#### 6. Explanation of symbols

 **0123** This product complies with Council Directive 93/42/EEC of 5 September 2007 concerning medical devices and bears the mark CE 0123 (TÜV SÜD Product Service GmbH).

Degree of protection against electric shock TYPE BF

Consult the instructions for use

Keep dry

Device serial number

Reference number = item number

**IP20**

Protected against solid foreign particles with a diameter of more than 12.5 mm, no protection against water.

Manufacturer

#### 7. Disposal

Batteries and technical appliances must not be disposed of with domestic waste, but should be handed in at the appropriate collection and disposal points.

- Keep the device clean. Check for cleanliness after use. Please use a soft, dry cloth for cleaning. Do not use benzene, thinners or other strong solvents.
- The cuff may absorb sweat and other fluids. Inspect the cuff for marks or discolourations after each use. Please use a synthetic detergent for cleaning and carefully wipe the surface. Do not scrub or machine-wash. Air-dry thoroughly.
- Make sure when storing the unit that no heavy objects are left lying on the unit.
- Do not apply force when folding or bending the cuff.
- If the unit has been stored at temperatures below 0 °C, let it stand at room temperature for at least 1 hour before use.
- If the device is not to be operated for an extended period of time, please remove the batteries. Store the batteries out of the reach of children.

The device has been manufactured and tested with all due care. Nevertheless, in the event of defects upon delivery, we provide a warranty under the following conditions:

During the warranty period of 3 years from the date of purchase, we remedy such defects at our discretion and at our expense in our plant through repair work or replacement delivery of a defect-free device.

The warranty does not cover normal wear and tear of wearing parts or damage caused by failure to observe the instructions for use, improper handling (e.g. unsuitable power sources, breakage, leaking batteries) and/or disassembly of the device by the buyer. Furthermore, the warranty does not constitute grounds for asserting claims for damages against us.

Warranty claims can only be asserted during the warranty period and upon presentation of the proof of purchase. In the case of a warranty claim, the device along with the proof of purchase and a description of the complaint must be sent to:

UEBE Medical GmbH  
Service-Center  
Zum Ottersberg 9  
97877 Wertheim, Germany

The cost of returning the device to our factory shall be borne by the sender. Complaints that are sent back without prepayment will not be accepted by UEBE.

The statutory claims and rights of the buyer against the seller (for example, defect claims, producer liability) are not restricted by this warranty.

**Please note: In the case of a warranty claim, please make sure to enclose the proof of purchase.**

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© Copyright 2018 UEBE Medical GmbH