ENG Please read the instructions carefully before use!

[INTENDED USE] Humasis COVID-19/Flu Ag Combo Test is one step in vitro diagnostic test based on an immunochromatographic assay. It is designed for qualitative detection of SARS-CoV-2 antigens and influenza A/B antigens in nasopharyngeal swab specimen of suspected patients.


[SUMMARY AND EXPLANATION] Coronavirus is a group of viruses that belongs to the Family Coronaviridae; a type of RNA virus of 27-32kb commonly found in birds and mammals including human. Coronavirus disease 2019 (COVID-19) is a disease caused by the severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) known to have originated from Wuhan city of China in December 2019, which has been declared a public health emergency.

[TEST PRINCIPLE] COVID-19 Ag Test A nitrocellulose membrane strip on the left side of the device contains one test line and one control line. The test line is pre-coated with anti-mouse monoclonal antibody to SARS-CoV-2 nucleocapsid and RBD for detection of SARS-CoV-2 antigens, and the control line is coated with goat anti-mouse IgG.


[CONTENTS] 25 Test devices packaged individually in aluminum pouch (25test/box) Disposable test tube with extraction buffer (25ea/box) Filter cap (25ea/box) Sterilized swabs for specimen collection (25ea/box) Instruction for use (1ea)

[STORAGE AND SHELF-LIFE] 18 months from manufacturing date at room temperature (2°C -30°C).


[TEST PROCEDURE] 1. Specimen collection: Use the swab included in the package to collect nasopharyngeal specimen. 2. Collected specimen should be tested immediately after collection for best result.

[INTERPRETATION OF RESULT] 1. Negative If no colored line appears in the test line (T) in both result windows and a colored line is present on the control region (C) in both result windows, then the result is negative.

3. Invalid The result is valid only when colored line is visible in the control line (C) in both result windows. If any of the control line (C) in both windows is absent, the test is invalid.

[PERFORMANCE CHARACTERISTICS] Limit of detection: The limit of detection (LoD) of Humasis COVID-19/Flu Ag Combo Test for COVID-19 virus is 5x10^2.7 TCID50/mL for influenza A H1N1 antigen is 2.075ng/mL, for influenza A H3N2 antigen is 5.5ng/mL, and for influenza B virus is 78ng/mL.

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Table listing various bacteria and viruses tested for cross-reactivity, including Mycoplasma pneumoniae, Streptococcus pneumoniae, Candida albicans, and others.

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Table listing interfering substances and their concentrations, such as Whole blood, Mucin, Chlorazepate, and others.

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Table showing clinical sensitivity and specificity data for COVID-19 Ag test in prospective study: Sensibilité clinique: 100.0% (5/5) (95% CI: 56.6%-100.0%) / Spécificité clinique: 100.0% (4/4) (95% CI: 91.4%-100.0%)

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[RÉFÉRENCES] 1) Zhu N, Zhang D, Wang W, et al. A Novel Coronavirus from Patients with Pneumonia in China, 2019. N Engl J Med 2020. 2) Huang C, Wang Y, Li X, et al. Clinical Features of Patients Infected with 2019 novel coronavirus in Wuhan, China. Lancet 2020. 3) Kang GK, Song KH, Cho PG, et al. Clinical and Epidemiological Characteristics of Spreaders of Middle East Respiratory Syndrome Coronavirus during the 2015 Outbreak in Korea. J Korean Med Sci 2017; 32:744-9.

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[REFERÊNCIAS] 1) Zhu N, Zhang D, Wang W, et al. A Novel Coronavirus from Patients with Pneumonia in China, 2019. N Engl J Med 2020. 2) Huang C, Wang Y, Li X, et al. Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet 2020. 3) Kang GK, Song KH, Cho PG, et al. Clinical and Epidemiological Characteristics of Spreaders of Middle East Respiratory Syndrome Coronavirus during the 2015 Outbreak in Korea. J Korean Med Sci 2017; 32:744-9.

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[LITERATURHINWEISE]

- [1] Zhu N, Zhang D, Wang W, et al. A Novel Coronavirus from Patients with Pneumonia in China, 2019. N Engl J Med 2020.
- [2] Huang C, Wang Y, Li X, et al. Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet 2020.
- [3] Kang DK, Song KH, Cho JG, et al. Clinical and Epidemiologic Characteristics of Spreaders of Middle East Respiratory Syndrome Coronavirus during the 2015 Outbreak in Korea. J Korean Med Soc 2017; 32:744-9.
- [4] WHO. Novel Coronavirus (2019-nCoV) situation reports. Verfügbar unter: https://www.who.int/emergencies/diseases/novel-coronavirus-2019/situationreports/ (Agerufen am 2. Feb. 2020).
- [5] Kim SH, Huh JH, Bae SY, Kim JS, Yoon SY, Lim CS, et al. The Dynamics of Respiratory Virus Infection in 2004-2006. The Journal of the Korean Society of Medical Diagnosis.

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