# Pe-Ha-Luron® F

Natriumhyaluronat / Sodium hyaluronate / Hialuronato sódico / Hyaluronate sodique

zum einmaligen Gebrauch  
for single use  
de un solo uso  
pour l'usage unique

Zusammensetzung: 1 ml enthält:	Content: 1 ml contains:	Composición: 1 ml contiene:	Composition : 1 ml contient	Pe-Ha-Luron®F 1.0 %	Pe-Ha-Luron®F 1.4 %	Pe-Ha-Luron®F 1.6 %	Pe-Ha-Luron®F 1.8 %	Pe-Ha-Luron®F 3.0 %
Natriumhyaluronat	Sodium hyaluronate	Hialuronato sódico	Hyaluronate sodique	10,000 mg	14,000 mg	16,000 mg	18,000 mg	30,000 mg
Natriumchlorid	Sodium chloride	Cloruro sódico	Chlorure de sodium	8,500 mg	8,500 mg	8,500 mg	8,500 mg	8,500 mg
Dinatriumhydrogenphosphat	Disodium hydrogenphosphate	Fosfato disódico	Hydrogénophosphate disodique	0,563 mg	0,563 mg	0,563 mg	0,563 mg	0,563 mg
Natriumhydrogenphosphat	Sodium dihydrogenphosphate	Fosfato sódico	Hydrogénophosphate sodique	0,045 mg	0,045 mg	0,045 mg	0,045 mg	0,045 mg
Wasser für Injektionszwecke	Water for injection	Agua para la inyección	Eau pour injection	q.s	q.s	q.s	q.s	q.s

Inhalt: Sterile Einmalspritze mit 1 ml 1,0% / 1,4% / 1,6% / 1,8% / 3,0% Natriumhyaluronat-Lösung und Einmal-Injektionskanüle, Luer-Lock  
 Content: Sterile single-use syringe with 1 ml 1,0% / 1,4% / 1,6% / 1,8% / 3,0% sodium hyaluronate solution and single-use injection cannula, Luer-Lock  
 Contenido: Jeringuilla estéril desechable de 1 ml Disolución al 1,0% / 1,4% / 1,6% / 1,8% / 3,0% de hialuronato sódico y cánula de inyección desechable, Luer-Lock  
 Contenu : 1 seringue stérile à usage unique de 1 ml, contenant une solution de 1,0% / 1,4% / 1,6% / 1,8% / 3,0% d'hyaluronate sodique, et une aiguille d'injection, Luer-Lock

## Gebrauchsinformation

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## Instruction for use

### Description:

Pe-Ha-Luron®F is a sterile, isotonic, pyrogen-free, viscoelastic preparation of highly purified sodium hyaluronate, dissolved in a buffered physiological solution (pH 6.8 to 7.5) for intraocular application. It is a clear solution in a ready-to-use syringe for single use. It is packaged in a sterile blister or edge-sealed pouch.

### Characteristics:

Sodium hyaluronate is a biopolymer obtained from bacterial biofermentation that belongs to the class of glycosaminoglycans. It has a high molecular weight, is highly purified and possesses favorable physical and rheological properties for ophthalmic surgery.

Hyaluronic acid is a biological and biocompatible substance.

The viscoelastic properties of Pe-Ha-Luron®F allow lubrication, support and protection of ocular tissues during ophthalmic surgery. It maintains the depth of the anterior chamber.

Pe-Ha-Luron®F forms a thin protective layer on the ocular cells and tissues and facilitates the insertion of the intraocular lens.

Pe-Ha-Luron®F is free from preservatives.

Due to its solubility in water, Pe-Ha-Luron®F can be easily removed by irrigation and aspiration at the end of surgery.

### Indications:

Pe-Ha-Luron®F serves as an adjuvant for the following operations:

- Anterior chamber surgery of the eye, and especially:
- Cataract surgery with or without intraocular lens implantation
- Glaucoma surgery

### Contraindications:

Hypersensitivity against sodium hyaluronate.

### Mode and duration of application:

If Pe-Ha-Luron®F was stored in the refrigerator, the syringe should be removed at least 20 minutes prior to use from the refrigerator. When the air has escaped from the ophthalmological Luer-Lock cannula, Pe-Ha-Luron®F can be carefully injected

## EN

into the anterior chamber using the enclosed sterile single-use cannula or a comparable cannula. Dosage depends on the kind of surgery. To replace loss of viscoelastic by outflow or irrigation, several injections can be made.

Remove Pe-Ha-Luron®F completely at the end of surgery using a suitable irrigation / aspiration device.

### Precautions:

Observe all usual precautions required during ophthalmic surgery!

Together with the enclosed sterile single-use cannula Pe-Ha-Luron®F is designed for single intraocular application only.

Check the integrity of the primary package (sterile pack) prior to the use of Pe-Ha-Luron®F to guarantee sterility.

After the intraoperative use of Pe-Ha-Luron®F a transient post-operative intraocular pressure increase may occur. A suitable therapy for reduction of the intraocular pressure decreases the pressure increase.

### Interactions with other devices or drugs:

There is an incompatibility of sodium hyaluronate and quaternary ammonium compounds such as solutions containing benzalkonium chloride. Therefore Pe-Ha-Luron®F must never get into contact with surgical instruments rinsed with these solutions or with ophthalmics containing quaternary ammonium compounds as preservatives.

### Shelflife:

Do not use Pe-Ha-Luron®F after the expiration date.

### Storage:

Store Pe-Ha-Luron®F between 2° and 25°C protected from light.

Do not freeze Pe-Ha-Luron®F - Avoid shock.

### Note:

For use during surgery:

Do not exert any more pressure on the plunger prior to removing the cannula from the eye. This is to avoid aspiration of air bubbles into the cannula.

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ALBOMED GmbH