NL

Belangrijk

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EN

Important

Please read this user’s manual carefully before using the device, and keep it for future reference. The charger is not included in the scope of supply. The medicsBlu® OraLightStick® is patent pending. The medicsBlu® OraLightStick® is not classified according the medical products law.

Danger

- Keep the charger away from water.
- Do not place it above or directly next to a bathtub or sink filled with water.
- Do not immerse the device in water or other liquids.
- After cleaning, make sure that the charger is completely dry before reconnecting it to the mains.

Warning

- Do not look directly at the blue light for a long period of time!
- The battery charging time is approx. 4 hours
- If it heats up too much during use, please switch it to the pulsing setting by pressing the on-switch again.
- Do not charge the device outdoors or near hot surfaces.
- Do not use the device if it is damaged (device or charging cable). This device does not contain any parts that the user can repair. If the device is damaged, please contact RP Medics (see title page International Warranty (warranty conditions)).
- This device can be used by children and persons with reduced physical, sensory or mental abilities or a lack of experience and knowledge provided they are supervised during use or have been instructed on how to use the device safely and understand the associated hazards. Cleaning and maintenance of the device must not be carried out by children without supervision.
- Children must not play with the device.

Attention

- Do not clean the device or the silicone hood in the dishwasher.
- After each use, disinfect the silicone hood of the medicsBlu® OraLightStick® for 5 minutes in boiling water.
- Replace the silicone hood every 6 months, or sooner if there are signs of wear and tear.
- Never use the medicsBlu® OraLightStick® if the silicone hood is damaged.
- The medicsBlu® OraLightStick® meets the safety standards for electromagnetic devices. If you have a pacemaker or another type of implant, consult your doctor or the manufacturer of the implanted device before using this device.

- If you have concerns about your health, consult your doctor before using the medicsBlu® OraLightStick® .
- This device is intended exclusively for cleaning the mouth and throat area. Do not use it for other purposes. If you experience discomfort or pain, stop using the device and consult your doctor.
- If you have had surgery on your teeth or gums in the past 2 months, you should consult your dentist before using the medicsBlu® OraLightStick®.
- The medicsBlu® OraLightStick® is a device for individual oral care and is not suitable for use by multiple patients in a medical practice or clinic.

Storage

If you decide not to use the medicsBlu® OraLightStick® for a long period of time, disconnect the USB charging cable, clean the medicsBlu® OraLightStick® and store it in a cool, dry place that is not exposed to direct sunlight.

Recycling

This product must not be disposed of with your regular household waste (2012/19/EU).

The product contains a battery that must not be disposed of with your regular household waste (2006/66/EC). Please hand the product in at an official collection point so that the battery can be removed professionally.

Used devices/used rechargeable batteries/normal batteries can be handed in free of charge at suitable collection points.

They are collected and sorted in the correct manner at these collection points, or prepared for reuse. Used devices/used rechargeable batteries/normal batteries may contain pollutants that are hazardous to the environment and human health. By recycling the raw materials contained in used devices/used rechargeable batteries/normal batteries, it is possible to make a contribution to environmental protection. Information for consumers in Germany: The return process detailed in point 1 is required by law. Used devices containing used batteries and/or used rechargeable batteries that are not enclosed by the used device must be detached from the battery before being handed in at a collection point. This does not apply if the used devices are still suitable for re-use (i.e. fully functional) and are separated from the other used devices. Collection and return points in Germany: https://www.stiftung-ear.de/; batteries can also be returned to retailers.

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RP Medics GmbH